
10055791

1/15






FIG. 1B





74




FIG. 2

_OVED	O.G. F	IG.
BY	CLASS	SUBCL
RAFTSMAN		

System	ECG	SpO2	NIBP	TEMP
Controls for Current Pr	or System S ofile	Setup		Monitor
	Fact	ory		Status
Alarm Sus	pend Time	180	▼seconds	Profiles
Alarm Vol	ume	High	ᢦ	
Tone Volu	me	6	ᢦ	
Heart Rate	e Sourse	AUTO [⊽	
				Color
				Cancel
				Save

System Setup Control Popup

FIG. 3A



	O.G. F	
Β̈́Υ	CLASS	SUBCL
RAFTSMAN		

System	ECG	SpO2	NIBP	TEMP
Controls for Current P	or System S rofile	Setup		Monitor
	Fact	ory		Status
Alarm Sus	spend Time	180 [seconds	Profiles
Alarm Vol	ume	180/// 150	4	
Tone Volu	ıme	120		
Heart Rate	e Source	90 60	i	
		30	ĺ	
	·		_	Color
				Cancel
				Save


Alarm Suspend Pull Down Has Been Clicked On

FIG. 3B



FIG. 4A

Screen format Control Panel



User Has Clicked on ECG



FIG. 4B

LOOSI791.OSOSO

Controls for	Alarm Se	 etup			
	Current		Low		
Parameter	Value	_Limit	Limit	Priority	
ECG	78	110	40	High	
SpO2	98	OFF	85	High	
NIBP Sys Dias	130 85	190 110	OFF OFF	High	
Mean	(110)	OFF	OFF		
TEMP	37.0	OFF	OFF	Medium	
☑ Display L	.imits?				Close

Alarm Setup Control Panel

FIG. 5



System ECG	SpO2	NIBP	TEMP
Controls for ECG			
Alarm High	110	$\overline{\mathbf{v}}$	
Alarm Low	40 [$\overline{\mathbf{v}}$	·
Alarm Priority	High [₹	
Size	1	▼ mV/cm	
Lead Select	11	$\overline{\Box}$	
Filter	60Hz	$\overline{\Box}$	• .
Sweep Speed	25 [▼ mm/sec	Color
☐Patient Paced			C0101
☑QRS Tone	Record	l On Alarm	Cancel
☐Display Units			
☑Cascade ECG			Save

ECG Popup Control Panel

FIG. 6A



OVED	O.G. F	IG.
BY	CLASS	SUBC'_
RAFTSMAN		




FIG. 6B




System ECG	SpC)2	NIBP	TEMP
Controls for Pulse	Oximetr	y		
Alarm High	OFF	♥		
Alarm Low	85	P		
Alarm Priority	High	P		
Average over	2	Þ	seconds	
Sweep Speed	25	₹	mm/sec	
☑ Display Pulse R				Calar
Display Pulse B	aı			Color
☐Display Units				Cancel
				Save

Control Popup for SpO₂

FIG. 7





System	ECG	SpC	2	NIBP] TEMP	
Controls for Alarm Limit		asive Blo	ood Pi Lov			
Systolic	19	0 🗖	OF	F ▽		
Diastolic	110	0 ~	OF	F ₹		
Mean	OF	F	OF	F ₹		
Alarm Prio	rity Hig	ıh ▽		ſ		
Displayed Values	SD	M &		ļ	Stat NIBPs for 5	
Auto Mode	OF	F $lacksquare$		į	Minutes	
1st Inflation	า 160) 🔻			Color	
☑Display	Units			ĺ		
	ion Tone			Į	Cancel	
					Save	

Control Popup for NIBP

FIG. 8



OVED O.G. FIG.

System	ECG	SpO2	NIBP	TEMP
Controls f	or Temper	rature		•
Alarm Hig	h L	OFF		
Alarm Lov	v [OFF		
Alarm Pri	ority [Medium		
Units		Celsius 🔻		
□Display	Units			
				·
				Color
			i	
				Cancel
				Save
	·			


Temperature Control Panel Popup

FIG. 9



1.0

ECG Processing Flow



G - ECG HR is 0(possible asystole or vfib) H - Beat Detection Signal with QRS information F - Filtred ECG Data for Beat Detection A - Raw ECG Data from A/D Converter B - Lead Selection Information C - Lead Selection Control Lines D - Valid ECG Data E - ECG Waveform Display Data

I - ECG HR Value

AUG 0 5 2002 REAL MAIN

CLOS SULLE

SMS Patient Data	a Transfer X
Temperature °F	98.6 Resp Rate (0-50)rpm
Pulse Rate(bpm)	78 Height (0-84)in 🔽
O2 Sat (%)	98 Weight (0-300)lbs
Systolic BP (mmHg)	130 Pain Scale (0-10)
Diastolic BP	Violence Scale (0-60)
	Store Patient Vitals Refresh Values Cancel

FIG. 12



RAFTSMAN

15/15

⊞SMS integrity Monit	HSMS integrity Monitor v2.0.0.1 (Simulation Mode)				_ EX
	Alarms Off			LEADII	<u>S</u>
78					Screen
Eag			-	1m/ C	Setup
Đ] } } }		}]		
ECG					Alarm Setup
	Alarms Off		<		
98 Sp02					Freeze
Pulse Rate=108		<i>}</i>	F	} }	Alarm Suspend
130	Time Rate NIBR 10:37a 118130/85/(110)	98.6∘F			
NIBP 85	80 128/81/(106) 84 125/81/(108) 77 122/78(105)	TEMP			Start NIBP
	10:37a 118 130/85(110)				Start
mmHg Alarms Off		Alarms Off	10:38 TIMER	10:38 am Dec 5, 2001 MER	- I